



**Cotton-top Tamarin** (*Saguinus oedipus*)

**Size:** Head and body length, 22.5 - 24 cm; tail, 33 - 40 cm **Weight:** 411 - 430 g

**Habitat:** Endemic to northwestern Colombia; found in rainforest, deciduous forest and secondary forest

**Surviving number:** Estimated at 6,000



Photographed by Lisa Hoffner

# WILDLIFE AS CANON SEES IT

Share and share alike. It's a motto cotton-top tamarins embody, as they share much more than a striking white topknot and mane. They share food liberally, especially fathers with offspring and older siblings with younger ones. They divvy up guard duty, with one keeping watch while the others forage. And they work together in defending their turf with vocal and physical displays; males fluff up their hair menacingly, approach invaders

and attack, while females scent-mark territory. But these cooperative creatures also share the perils of habitat loss and capture for the pet trade, which are taking a heavy toll on populations.

As we see it, we can help make the world a better place. Raising awareness of endangered species is just one of the ways we at Canon are taking action—for the good of the planet we call home. Visit [ngm.com/canonwildlife](http://ngm.com/canonwildlife) to learn more.