



White-winged Guan (*Penelope albipennis*)

Size: Head and body length, approx. 50 cm; tail, approx. 35 cm **Weight:** Approx. 1.6 kg **Habitat:** Dry wooded slopes and ravines on the west slope of the Andes in Peru **Surviving number:** Estimated at 200



Photographed by Roland Seitre

WILDLIFE AS CANON SEES IT

Now hear this. Announcing his claim to a territory, the male white-winged guan calls with a distinctive “jar jar jar” sound, gradually increasing in speed and ending in a quick “ha ha ha.” Monogamous couples share the duty of tending the nest, and young are able to venture out soon after hatching to forage for leaves, buds, fruits, flowers and seeds. Even so, their parents continue to help feed and care for them until they are a year old. But it

is not the young alone that need help to survive; once feared extinct and still critically endangered, the species is threatened by habitat loss and the persistence of opportunistic poaching.

As we see it, we can help make the world a better place. Raising awareness of endangered species is just one of the ways we at Canon are taking action—for the good of the planet we call home. Visit canon.com/environment to learn more.

Canon